



IT'S SIZE IS NO
MATCH FOR ITS POWER

iStir Uno+

(Motorized Ultra Thin Magnetic Stirrer)

iStir Uno+

(Motorized Ultra
Thin Magnetic Stirrer)

Immense power in a compact package, iStir Uno+ promises reliable performance that is outclassed only by its durability.

Made from the best quality materials iStir Uno+ makes working in the lab much easier.

With convenience built right into the design and user-interface of the device, working with iStir Uno+ is an amazing experience.

Make the most of your time with timesaving features.

This a device that is sure to offer consistent performance that will last for years.



More Than Meets The Eye

Meet, iStir Uno+, a class-leading motorized magnetic stirrer that promises exceptional performance along with unmatched durability. The device's entire production cycle including design, manufacturing and quality control is extremely stringent to ensure that the end-user's experience with the product is satisfactory. As a result, the device is full of features that will make the user's work easier by asking for less effort and time.



BLDC Motor



Stainless Steel Top Surface



Microprocessor Controlled

Pulse Mode



Gradual Acceleration



IP44 Compliant Design



Experience Brilliance!

iStir Uno+ is a brilliant device that is made to ensure that you a superlative experience in the lab. Designed to meet all of your stirring needs, the device incorporates brilliant performance in a small package all the while delivering exceptional ease of use.



Last Run Memory Function



Adjustable Speed Setting

Variable Timer Setting



Low Aspect Ratio



Specifications

Specifications	iStir Uno+
Max. Capacity	10 litre
Motor type	BLDC
Speed setting	200-2200 RPM (every 100 RPM)
Run time	1-999 mins & infinite mode
Recommended stir bar length	35 mm
Display	Digital-shows real time speed
Pulse mode	Yes
Last run memory	yes
Dimensions (W X D X H)	335 x 200 x 29
Protection class	IP 44 compliant design (water spill resistant)
Weight	2.5 kg
Standard accessory	1 high strength stir bar
Protection class	IP 44 (water spill resistant)

Ordering Information

iStir Uno+		
Product Name	Cat #	Details
iStir UNO+	MS02-1110 UK	Single position UNO+, all SS stirrer with UK pin
iStir UNO+	MS02-B4	Stir bar 35 mm length

NEUATION
TECHNOLOGIES
 www.neuation.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.